



# *Precision-Fit*

Dentures



**EVERY DENTIST'S FIRST GOAL IS TO PROTECT THE HEALTH OF YOUR NATURAL TEETH.**

Yet, circumstances beyond your control sometimes lead to having to face the reality of complete dentures. Dr. Clay Couvillon of Sullivan Dental continually educates himself and his team in the latest denture techniques and strategies to make every patient's experience a positive and successful one.



The Couvillon team provides confident smiles that can make patients look younger, as well as give them the ability to eat healthy foods. The team works hard to enhance the quality of life of each patient that walks through its doors. Life is too short to have your teeth cause you embarrassment or discomfort. Now is the time to enjoy your smile!

"I dedicate myself to each patient's satisfaction. We customize a process to ensure that your new dentures will fit comfortably and increase your self confidence and quality of life."

---

*Dr. Clay Couvillon, D.D.S.*

## What Makes Couvillon Precision-Fit Dentures Unique

We gave our dentures the Couvillon Precision-Fit name for their handcrafted construction which suits each patient's individual appearance, needs, and precise function.

### Customized Fit

Customizing your dentures begins with a two-step impression process. This process gives us an accurate outline of the surrounding tissue on the floor of the mouth and cheeks. For you, these impressions translate to a better, more comfortable fit.

### Precise Function

You will want your new dentures to stay in place and allow you to eat foods you thought you could no longer enjoy. We focus next on aligning and balancing the bite using two state-of-the-art techniques.



Macia André's before and after photographs. You can see the natural "facelift" effect achieved with these dentures.

### A Natural "Facelift"


We design your dentures to fit the cosmetic proportions of your unique anatomy. Dr. Couvillon's precise measurements help bring back your genuine smile, providing a natural "facelift" and restoring a more youthful look.

### Dentures for Life

Creating quality dentures takes time and a dedication to the final product. Our Precision-Fit process includes a duplicate denture that will allow future remakes to be less costly and time-consuming for you.

"I wasn't sure what to expect when I walked through the doors of Sullivan Dental Center. During my consult with Dr. Couvillon and his team, they truly listened to my concerns and made it easy to decide on an implant-retained denture. He treated me like family, and if something wasn't right, he and his team worked until it was perfect. Today, I have a beautiful smile I can be proud of."

*Ms. Macia André*





## Customized by Couvillon for Comfort and Function

### What makes Couvillon Precision-Fit Dentures different from the rest?

- Dr. Couvillon is the **first North American dentist trained in mandibular suction**, a technique he uses to take impressions that ensure your best fit and maximum comfort.
- To maximize the function of your new smile, we customize your dentures using techniques that **precisely align and balance** the upper and lower teeth.
- Confidence in your appearance means a lot to you, and to Dr. Couvillon. His precision crafting leads to a more youthful face proportion, **like a facelift without the surgery**.
- **We take the time to get it right.** Each Precision-Fit Denture has duplication capabilities so that future remakes are less costly and time-consuming for you.

## Dr. Couvillon: One of a Kind

The first general dentist in North America to train in mandibular suction, Dr. Couvillon learned this revolutionary technique from Japanese specialist Dr. Jiro Abe in 2016. Dr. Couvillon uses mandibular suction to improve the fit of all types of dentures on all patients. He finds the technique particularly helpful with securely fitting lower dentures without the aid of implants.



“I am very pleased with Dr. Couvillon and his team. I would highly recommend him to anyone who is thinking about getting new dentures.”

*Mr Douglas Rodney*



## Couvillon Precision-Fit Dentures: Three Options

### Tissue-Supported

Tissue-supported dentures are removable dental devices to replace your natural teeth. The upper and lower dentures fit onto the gums, secured by suction. Dr. Couvillon's improved fitting techniques make this the right choice for patients who are not candidates for dental implants.



### Implant-Supported

Also removable, implant-supported dentures replace your natural teeth by anchoring to dental implants inserted into the jawbone. Wearers remove the dentures to clean them. The implants give the dentures structural support and stability, and dramatically reduce bone loss over time. This option is the most cost-effective way to give your new teeth a more natural feel and open the palate to improve taste.

### Implant-Retained

Implant-retained dentures differ from implant-supported dentures because they are fixed dental devices. Wearers maintain the dentures like their natural teeth without removing them. They have the same structural support and stability and bone preservation benefits. No one else has to know you're wearing dentures when you choose the implant-retained option.



## What to Expect Free Initial Consult and Oral Exam

**Dr. Couvillon will meet personally** with you to talk about your needs and goals. He will complete a free oral examination to evaluate your dental bone and tissue.

## CBCT 3D Imaging Scan

**Dr. Couvillon's team** uses three-dimensional imaging to help him get a more accurate view. This will help determine the most beneficial dentures for you.

## Developing Your Treatment Plan

**Once you decide** on the best type of dentures for your needs, one of our specialized treatment coordinators will personally assist you throughout your entire experience.



## The Denture Appointment Process

**All denture patients visit for two** appointments that include our two-step impression process, followed by a minimum of two wax try-in appointments. We want to be sure you love the fit and function of your final dentures.

For implant-supported and implant-retained denture patients, the process begins instead with two appointments for initial impressions and placement of implants. These patients require about 4-6 months of healing time before denture fittings can begin.

**All of Our Precision-Fit Dentures Require a Minium of 4 Appointments**



FEATURES	TISSUE-SUPPORTED DENTURE	IMPLANT-SUPPORTED DENTURE	IMPLANT-RETAINED DENTURE
NATURAL SMILE	✓✓✓	✓✓✓	✓✓✓
RESTS ON GUMS	✓✓✓	✓✓	--
NORMAL PRESSURE ON GUMS	✓✓✓	✓✓	--
BITE WITH MORE FORCE	✓✓	✓✓✓	✓✓✓
REMOVE TO CLEAN	✓✓✓	✓✓	✓
REQUIRE PERIODIC ADJUSTMENTS, MAINTENANCE	✓✓✓	--	--

✓✓✓ = YES    ✓✓ VERY LITTLE    -- = NONE

**Life is too short** to have your teeth cause you embarrassment or discomfort. Now is the time to enjoy your smile! Call us today at **225-635-4422** or send us a message at [clay@drclaycouvillon.com](mailto:clay@drclaycouvillon.com) to schedule a free consultation and exam appointment.

Now Is the Time  
to Enjoy Your Smile!

# COUVILLON

*Dentistry Redefined*



5436 Commerce Street  
St. Francisville, LA 70775

225.635.4422 office  
225-635-2171 fax

[sullivandental.com](http://sullivandental.com)  
[sullivandental@bellsouth.net](mailto:sullivandental@bellsouth.net)